

DSB Analysis in rDNA

1. Suspend 5×10^7 cells in 50 μ l of cold L buffer (0.1 M EDTA, pH 8, 0.01 M Tris, pH 7.6, and 0.02 M NaCl) in an Eppendorf tube containing 5 μ l of cold zymolyase (zymolyase-20T, *Arthrobacter luteus*, 20000 U/gm; ICN) at a stock concentration of 20 mg/ml.
2. Warm cell suspension by incubating the tube in a 42° C water bath for a few seconds and mix briefly with 50 μ l of 1% low melting agarose prepared in L buffer cooled down to 42° C.
3. Placed mixture in an ice water bath.
4. Transfer gel plug to a 20 ml round-bottomed glass tube containing 5 ml of 0.5 M EDTA, pH 8.0 and 0.01 M Tris, pH 7.6 plus 1% β -mercaptoethanol. (The gel plug can be easily flushed out into the glass tube with the buffer).
5. Incubate the gel plug at 37° C for 24 h.
6. Remove buffer and add 0.5 ml of the same buffer containing 0.5 mg/ml zymolyase.
7. Incubate the gel plug at 37° C for 24 h.
8. Remove buffer and add 1 ml of L buffer containing 1 mg/ml proteinase K and 2.5% Sarkosyl.
9. Incubate the gel plug at 50° C for 24 h.
10. Cool gel plug on ice, remove buffer and add 0.5 ml of the same buffer containing proteinase K and Sarkosyl.
11. Incubate at 50° C for 24 h.
12. Incubate the gel plug in 5 ml of Tris-EDTA (pH 7.6) containing 40 μ g/ml phenylmethylsulfonyl fluoride at 50° C for 1 h.
13. Incubate in 5 ml of Tris-EDTA at 50° C for 1 h.
14. Incubate the plug with 100 μ g/ml RNase A in 1 ml restriction enzyme buffer at 50° C for 3 h.
15. Remove buffer and incubate the plug in 0.5 ml of restriction enzyme buffer containing BglII (100 units) at 37° C overnight.
16. Soak the gel plug in 5 ml of Tris-EDTA at 4° C for 0.5 h.
17. Drain buffer, add 20 μ l of 6X loading buffer (28) to the gel plug and incubate on ice for 0.5 h.
18. Load the gel plug into the well of a 1% agarose gel and run in TBE buffer at 1.5-2 V/cm at 4° C for 24 h.
19. Stain gel in 1 μ g/ml ethidium bromide for 2-4 h, inspect under UV and use for Southern blot.

Probes and labelling-probe DNA was prepared by PCR from a rDNA plasmid with the following pairs of primers: probe A, ER5A-up, 5'-GCC ATT TAC AAA AAC ATA ACG-3' and ER5A-lower, 5'-GGG CCT AGT TTA GAG AGA AGT-3'; probe B, X35S-up, 5'-ATA TCA ACC CTG ACG GTA GAG-3' and X35S-lower, 5'-CAT GGT ATA ACT GTG GTA ATT CTA GAG-3'; probe C, BglII end-up, 5'-ACA GAT GTG CCG CCC CAG CCA AAC TCC-3' and BglII end-lower, 5'-CCT GGA TAT GGA TTC TTC ACG GTA ACG-3'. The PCR products were gel-purified. Random-primed labelling kit (Roche Applied Science GmbH) was used to label the probe DNA in the presence of four [α -³²P]dNTPs.